**Best Fit**

**08**

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SISTEMA CLICK & GRIP

**ISOGRIP®**

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/ 8 EASY STEPS

CLICK & GRIP SYSTEM  
FOR IMMEDIATE LOADING  
ABOUT PILAR MULTI-UNIT

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L

#### / OUTER CYLINDER

Bonds the Isogrip® System to the slipped prosthesis.

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Código: G09160

#### / FIX SCREW

Responsible for joining the Isogrip® System to the Multi-Unit abutment.

Code: G09160

#### / Carrying out immediate loads with the reversible Isogrip® System.

Carrying out immediate loads on any type of dental over Multi-Unit pilar, thanks to the Isogrip®.

This system allows:

- Make the screw channels in their exact position;
- Time saving;

In just 8 steps, we obtain a denture screwed on implants.

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#### / ACTIVATION SCREW

Piece that houses the seat of the fixing screw; Activates the ultimate expanding core when tighten with the Activation Key (F).

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#### / DEFINITIVE EXPANSION CORE

Piece that is attached to the provisional cylinder through an expansion of his upper area, produced by the pressure exerted by the activator screw.

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#### / PROVISIONAL EXPANSION CORE

Piece geometrically analogous to the Isogrip abutment, with the only difference that it is machined in PEEK and does not require an activator screw.

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## / THE ISOGRIP® SYSTEM FOR **IMMEDIATE LOAD.**

COMPLETE STRUCTURE CREATION.

### / HOW IT WORKS

#### /STEP 01 FIXATION

Screw the **Provisional Expander Core** (E) to the Multi-Unit pillars with the **Fix Screw** (B), and fit the **External Cylinders** (A) exerting slight pressure, until a “click” is felt.

#### /STEP 02 OVERFLOW

Overlay the prosthesis with resin.

#### /STEP 03 UNION

Union of the overlay prosthesis on the **Isogrip® Provisional System** (E) and solidification of the resin.

#### /STEP 04 SEPARATION

Separate the set -formed by the prosthesis including the **External Cylinders** (A) linked inside- and from the Multi-Unit pillars, exerting slight upward pressure.

#### /STEP 05 DRILLING

Fit the **Drilling Guides** (H) on the inside of the **External Cylinders** (A), and insert the **Cutting Drill** (G) to drill the channels of the **Fixation Screws** (B) in this exact position.

#### /STEP 06 REMOVAL

Unscrew the **Provisional Expander Nucleus** (E), from the Multi-Unit pillars.

#### /STEP 07 INSERTION

Insert one by one, the **Definitive Expansion Cores** (D) in the **External Cylinders** (A), and make sure that the **Activating Screw** (C) does not exert pressure on the upper part of the **Definitive Expansion Core** (D). Next, tighten the **Activation Screws** (C) with the **Activation Wrench** (F) in an counter clockwise.

#### /STEP 08 FINAL SCREWING


Position the **prosthesis** with the **Definitive Isogrip® System** included on the Multi-Unit pillars. Finish by screwing the complete structure using the **Fixing Screws** (B) to the implant.

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Código: E001564

**/ ACTIVATION KEY**

Key to tighten the **Activator Screw (C)**.

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Código: E000502

**/ CUTTING DRILL**

Drill required to drill the **Fix Screw (B) channels**.

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Código: E73001

**/ DRILLING GUIDE**

Guide to direct the cutter on this **exact position** and make the screw **channel super precisely**.

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Código: D30054

**/ SCREWDRIVER**

Used to screw the **Fix Screw (B)** to the **Isogrip® System**.

Code: D30054